

Polarization of light

Basic Information:




Polarization of EM wave




- Light is considered polarized along Electric field
- Unpolarized light has random polarization direction


An unpolarized light is polarized by a polaroid

The polarised wave may be blocked by making *pass-axis* of the polaroid perpendicular to direction of polarization





Malus' Law



- 1) The 1st polarizer is used to polarize unpolarised light in a plane
- 2) The 2^{nd} polarizer (analyzer) is rotated w.r.t. the 1^{st} polarizer by an angle θ

Unloparised light through polariser

The intensity of an unpolarized light across a plane polarizer also reduces following the relation $I_0 = I\cos^2\theta$, I is the intensity before polarizer


When averaged over all possible angles, the total intensity reduces by half

$$I_0 = I < \cos^2 \theta > = \frac{I}{2}$$


Polarization by reflection

Unpolarized incident light on the surface of a material

Dipoles oscillate with the E-field and emits radiation




Polarization by reflection: Brewster's law



$$\theta_B = tan^{-1} \left(\frac{n_2}{n_1} \right)$$

At this angle of incidence, a plane π -polarized light has zero reflection coefficient. So, for unpolarized incident ight, the reflected ray will be ane polarized and refracted ray will be partially polarized

Polarization by reflection: Brewster's law



$$n_1 \sin \theta_B = n_2 \sin \theta_r$$

$$n_1 \sin \theta_B = n_2 \sin \left(\frac{\pi}{2} - \theta_B\right) = n_2 \cos \theta_B$$

$$\tan \theta_B = \frac{n_1}{n_2}$$

Brewseter's angle

$$\theta_B = tan^{-1} \left(\frac{n_2}{n_1} \right)$$

At this angle of incident, a plane polarized light has zero reflection coefficient. So, for unpolarized incident light, the reflected ray will be plane polarized and refracted ray will be partially polarized

Examples of polarization by reflection





Look at the window!!



Superposition of two plane polarized wave

$$\vec{E} = \hat{i}E_x + \hat{j}E_y$$

$$E_x = E_{x0} \cos(kz - \omega t)$$

$$E_v = E_{v0} \cos(kz - \omega t + \delta)$$



$$\frac{E_{y}}{E_{y0}} = \cos(kz - \omega t)\cos \delta - \sin(kz - \omega t)\sin \delta$$

$$\frac{E_{y}}{E_{y0}} = \frac{E_{x}}{E_{x0}} \cos \delta - \sqrt{1 - \frac{E_{x}^{2}}{E_{x0}^{2}}} \sin \delta$$


$$\left(\frac{E_{y}}{E_{y0}} - \frac{E_{x}}{E_{x0}}\cos\delta\right)^{2} = \sin^{2}\delta - \frac{E_{x}^{2}}{E_{x0}^{2}}\sin^{2}\delta$$

$$\left[\frac{E_{y}^{2}}{E_{y0}^{2}} - 2\frac{E_{x}}{E_{x0}}\frac{E_{y}}{E_{y0}}\cos\delta + \frac{E_{x}^{2}}{E_{x0}^{2}} = \sin^{2}\delta\right]$$

$$\frac{E_{y}^{2}}{E_{y0}^{2}} - 2\frac{E_{x}}{E_{x0}}\frac{E_{y}}{E_{y0}}\cos\delta + \frac{E_{x}^{2}}{E_{x0}^{2}} = \sin^{2}\delta$$


This is an equation of ellipse whose major axis is making an angle say α

$$\tan 2\alpha = \frac{2E_{x0}E_{y0}\cos\delta}{E_{x0}^2 - E_{y0}^2}$$



Linearly polarized

$$\delta = 2m\pi$$
 $m = 0, 1, 2, 3....$



$$\frac{E_y}{E_x} = \frac{E_{y0}}{E_{x0}}$$

$$\theta = \tan^{-1} \left(\frac{E_{y0}}{E_{x0}} \right)$$


$$E_x = E_{x0} \cos(kz - \omega t)$$

$$E_{y} = E_{y0} \cos(kz - \omega t + \delta)$$

$$E_{y} = E_{y0} \cos(kz - \omega t + 2m\pi)$$

$$E_{y} = E_{y0} \cos(kz - \omega t)$$

$$\delta = (2m+1)\pi$$
 $m = 0,1,2,3...$



$$\frac{E_y}{E_x} = -\frac{E_{y0}}{E_{x0}}$$

$$\theta = -\tan^{-1}\left(\frac{E_{y0}}{E_{x0}}\right)$$

$$E_v = E_{v0} \cos(kz - \omega t + 2m\pi + \pi)$$

$$E_{y} = -E_{y0}\cos(kz - \omega t)$$

Elliptically polarized

$$\delta = (2m+1)\frac{\pi}{2}$$
 $m = 0,1,2,3...$ $\Rightarrow \frac{E_y^2}{E_{y0}^2} + \frac{E_x^2}{E_{x0}^2} = 1$



$$\frac{E_{y}^{2}}{E_{v0}^{2}} + \frac{E_{x}^{2}}{E_{x0}^{2}} = 1$$

At
$$z = 0$$


$$E_{x} = E_{x0} \cos(\omega t)$$

$$E_{v} = E_{v0} \cos(\omega t - \delta)$$

For $\delta = \frac{\pi}{2}$; $5\frac{\pi}{2}$; $9\frac{\pi}{2}$

$$E_x = E_{x0} \cos(\omega t)$$

$$E_{v} = E_{v0} \sin(\omega t)$$



 χ


Counter-clock wise rotation with time

 $E_{x0} = E_{v0}$


For
$$\delta = 3\frac{\pi}{2}$$
; $7\frac{\pi}{2}$; $11\frac{\pi}{2}$

$$E_x = E_{x0} \cos(\omega t)$$


$$E_{y} = -E_{y0}\sin(\omega t)$$




Clock wise rotation with time



$$E_{x0} = E_{y0}$$



A substance is <u>Optically Active</u> if it rotates the plane of polarized light.




Rotation is clockwise

Example: Glyceraldehyde

Rotation is counterclockwise

Example: D-fructose

Chiral molecules



Source of image: Wikipedia

There is no set of translation and rotations that can map the left-hand molecule into the right-hand side molecule

A solution containing one form of such asymmetric molecule may rotate the plane of polarized light

Specific rotation [S] of a chiral substance

If θ is the angle of rotation for a solution of concentration c and length L, where


- Length (*L*) measured in dm.
- θ is the angle of rotation.
- c measured in gm/cc

$$S = \frac{\theta}{Lc}$$

S is the rotation produced by a column of solution of length 1 decimeter and containing 1 gm of the active substance per cm³ of the solution at a particular **temperature and for a given wavelength of incident light:**


S is a function of temperature and wavelength

Half-Shade Plate





Quartz of thickness such that a path difference of $\frac{\lambda}{2}$ or a phase difference of π is introduced between two perpendicular components of electric field, one parallel to optic-axis and another perpendicular to it.


Other half made of glass, of adequate thickness so that light intensity is same as coming from the quartz side, when pass axis of analyzer (discussed in the slide 18) lies in the plane of polarization of light coming from half-shade plate.





Half-Shade Plate












Experimental Setup



Experiment to determine S

- Measure the change in analyzer pass-axis angle between distilled water and then with sugar solution
- Repeat it for several concentrations of sugar solution

Plot a graph of θ vs c.

<u>Table 3</u> Determination of the specific rotation of the solution from the plot			
(cm)	c (gm/cm ³)	θ (degree)	$s = \frac{10\theta}{lc}$ (degree cm ³ decimeter ⁻¹ gm ⁻¹)

S is determined